


CLIA

Product Name	Catalog #	# of Tests	Qualitative / Quantitative	Test Format	Sensitivity
--------------	-----------	------------	----------------------------	-------------	-------------

AccuDiag™ Cancer (Tumor Marker) Assays



AFP	9040-16	96	Quantitative	~65 min	2.0 ng/mL
CEA	9041-16	96	Quantitative	~65 min	1.0 ng/mL
PSA	9042-16	96	Quantitative	~125 min	0.5 ng/mL
Free Beta hCG	9043-16	96	Quantitative	~65 min	0.1 mIU/mL
Free PSA	9054-16	96	Quantitative	~125 min	0.05 ng/mL
Beta 2- Microglobulin	9076-16	96	Quantitative	~65 min	2.0 µg/mL
NSE	9056-16	96	Quantitative	~65 min	2.0 ng/mL
CA-12-5	9057-16	96	Quantitative	~125 min	5 U/mL
CA-19-9	9059-16	96	Quantitative	~125 min	5 U/mL
CA-15-3	9075-16	96	Quantitative	~125 min	5 U/mL
Ferritin	9019-16	96	Quantitative	~75 min	5.0 ng/mL

AccuDiag™ Cardio-Vascular Monitoring Assays



CRP	9053-16	96	Quantitative	~125 min	5 U/mL
-----	---------	----	--------------	----------	--------

CLIA

Product Name	Catalog #	# of Tests	Qualitative / Quantitative	Test Format	Sensitivity
--------------	-----------	------------	----------------------------	-------------	-------------

AccuDiag™ Allergy Assays



IgE	9050-16	96	Quantitative	~65 min	5 IU/mL
-----	---------	----	--------------	---------	---------

AccuDiag™ Fertility Assays



LH	9012-16	96	Quantitative	~65 min	2 mIU/mL
FSH	9013-16	96	Quantitative	~65 min	2.5mIU/mL
Prolactin (PRL)	9014-16	96	Quantitative	~65 min	2 ng/mL
hCG	9015-16	96	Quantitative	~65 min	2 mIU/mL

AccuDiag™ Steroids Assays

Progesterone	9025-16	96	Quantitative	~95 min	0.2 ng/mL
Estradiol	9028-16	96	Quantitative	~95 min	5 pg/mL
Testosterone	9052-16	96	Quantitative	~95 min	0.05 ng/mL

AccuDiag™ Thyroid Assays

T3	9001-16	96	Quantitative	~65 min	0.03 ng/mL
T4	9003-16	96	Quantitative	~50 min	0.5 µg/mL
TSH	9005-16	96	Quantitative	~65 min	0.03 ng/mL
Anti-Tg Min. Order: 5 kits	9010-11	96	Qualitative	~65 min	100%

CLIA

Product Name	Catalog #	# of Tests	Qualitative / Quantitative	Test Format	Sensitivity
AccuDiag™ Infectious Disease Assays					
H. pylori IgA Min. Order: 5 kits	9060-11	96	Qualitative	~65 min	N/A
H. pylori IgG Min. Order: 5 kits	9061-11	96	Qualitative	~65 min	99%
H. pylori IgM Min. Order: 5 kits	9062-11	96	Qualitative	~65 min	100%
Toxoplasma IgG Min. Order: 5 kits	9071-11	96	Qualitative	~65 min	94%
Toxoplasma IgM Min. Order: 5 kits	9072-11	96	Qualitative	~65 min	100%
CMV IgG Min. Order: 5 kits	9063-11	96	Quantitative	~65 min	97%
CMV IgM Min. Order: 5 kits	9064-11	96	Qualitative	~65 min	100%
HSV1 IgG Min. Order: 5 kits	9065-11	96	Qualitative	~65 min	100%
HSV1 IgM Min. Order: 5 kits	9066-11	96	Qualitative	~65 min	100%
HSV2 IgG Min. Order: 5 kits	9067-11	96	Qualitative	~65 min	100%
HSV2 IgM Min. Order: 5 kits	9068-11	96	Qualitative	~65 min	70%
Rubella IgG Min. Order: 5 kits	9069-11	96	Qualitative	~65 min	100%
Rubella IgM Min. Order: 5 kits	9070-11	96	Qualitative	~65 min	100%
HIV 1/2	9073-32	96	Qualitative	~65 min	N/A
Syphilis IgG/IgM	9074-32	96	Qualitative	~65 min	N/A
HCV	9075-32	96	Qualitative	~65 min	N/A
HBsAg Min Qty: 50 kits	9076-32	96	Semi- Quantitative	~65 min	N/A
HBsAb Min Qty: 50 kits	9077-32	96	Semi- Quantitative	~65 min	N/A
HBeAg Min Qty: 50 kits	9078-32	96	Semi- Quantitative	~65 min	N/A
HBeAb Min Qty: 50 kits	9079-32	96	Semi- Quantitative	~65 min	N/A
HBcAb Min Qty: 50 kits	9080-32	96	Semi- Quantitative	~65 min	N/A
AccuDiag™ Diabetes					
Insulin	9021-16	96	Quantitative	~65 min	2μIU/mL
C-peptide	9022-16	96	Quantitative	~65 min	0.05 ng/ml.

